



Member of LEONI Group



Autoclavable Fibers for Diode, Nd:YAG and Holmium Lasers

Type	Product Code	Core	Outer Diameter	Length	Color
200 µm	FTIR272/300ST-5/SL-F	272 µm	550 µm	5 meters	blue
365 µm	FTIR365/400ST-3/SL-F	365 µm	700 µm	3 meters	blue
550 µm	FTIR550/605ST-3/SL-F	550 µm	850 µm	3 meters	blue
1000 µm	FTIR1000/1100ST-3/SL-F	1000 µm	1450 µm	3 meters	blue

All our medical products are marked with CE 0197.

- Advantages:
- Products applicable up to 142 °C / 287,6 °F
 - applicable for Holmium Laser
 - high power F-SMA905 connector
 - accessories (stripper, cleaver) available

Further colors, diameters and lengths are available.

FiberTech produces special silica fibers with core diameters from 10 µm to 2000 µm with different numerical apertures, coatings and jackets. Fiber assemblies include cables for laser beam transmission, fiber tapers and fiber bundles. Please contact us for more information and for sample requests:

EN ISO 9001 + 13485 certified

FDA Reg. No. 3005128587

FIBERTECH® GmbH
 Nalepastraße 170 –171
 12459 Berlin
 Germany
 phone: [+49] 030 53 00 58-0
 fax: [+49] 030 53 00 58-58
 e-mail: medical@fibertech.de
 internet: www.fibertech.de
www.fibertechgroup.com

