




FIBERTECH®

SPECIAL OPTICAL  
FIBER TECHNOLOGIES

Member of LEONI Group



## Laserfiber for intravascular laser treatment

Type	Product Code	Core	Outer Diameter	Length	Color
400 µm	FTIR400T-3/SM-F-L	400 µm	750 µm	3 meters	transparent
600 µm	FTIR600T-3/SM-F-L	600 µm	980 µm	3 meters	transparent

All our medical products are marked with CE 0197.

- Advantages:
- sterile EO, double packed
  - IR-Fiber, low OH
  - SMA905 connector, cost-effective
  - accessories (stripper, cleaver) available

Further colors, diameters and lengths are available.

FiberTech produces special silica fibers with core diameters from 10 µm to 2000 µm with different numerical apertures, coatings and jackets. Fiber assemblies include cables for laser beam transmission, fiber tapers and fiber bundles. Please contact us for more information and for sample requests:

EN ISO 9001 + 13485 certified

FDA Reg. No. 3005128587

FIBERTECH® GmbH  
Nalepastraße 170 –171  
12459 Berlin  
Germany  
phone: [+49] 030 53 00 58-0  
fax: [+49] 030 53 00 58-58  
e-mail: [medical@fibertech.de](mailto:medical@fibertech.de)  
internet: [www.fibertech.de](http://www.fibertech.de)  
[www.fibertechgroup.com](http://www.fibertechgroup.com)

