



Ophthalmic Laser Probes

Type	Product Code	Core	Length	Color
Endoprobe, straight	FTOPH200-2,5/SES	200 µm	2,5 meters	white/blue
Endoprobe, curved illuminated	FTOPH200-2,5/SIC	200 µm	2,5 meters	white/blue
Cycloprobe, straight 900 µm Ball	FTOPH200-2,5/SCS	600 µm	2,5 meters	white/blue
Retinopexiprobe	FTOPH600-2,5/SES	600 µm	2,5 meters	white/blue

All our medical products are marked with CE 0197.

- Advantages:
- excellent mechanical and optical quality
 - sterile EO, double packed
 - IR-fiber, low OH
 - 20, 23, 25 Gauge available
 - angled probe with different bend radius available

Further colors, diameters and lengths are available.

FiberTech produces special silica fibers with core diameters from 10 µm to 2000 µm with different numerical apertures, coatings and jackets. Fiber assemblies include cables for laser beam transmission, fiber tapers and fiber bundles. Please contact us for more information and for sample requests:

EN ISO 9001 + 13485 certified

FDA Reg. No. 3005128587

FIBERTECH® GmbH
 Nalepastraße 170 –171
 12459 Berlin
 Germany
 phone: [+49] 030 53 00 58-0
 fax: [+49] 030 53 00 58-58
 e-mail: medical@fibertech.de
 internet: www.fibertech.de
www.fibertechgroup.com

